
# Documentazione guaina Tube assembly documentation Documentation gaine équipée

## C341 V

Nr. di matricola  
Tube No.  
Nr de série

CE 0051

Questa documentazione deve essere fornita all'utilizzatore del complesso tubo guaina  
The contents of this documentation must be transmitted to the user of the tube-assembly  
Le contenu de cette documentation doit être transmis à l'utilisateur de la gaine équipée

Documentazione N° Documentation N° N° de Documentation	Revisione Edition Version	Data di edizione Date of release Date de l'édition	Testo originale Original text Texte original
C341V	A	22.02.2022.	italiano / italian / italien

---

**Sommario - Table of contents - Table des matières**

---

Sommario - Table of contents - Table des matières .....	2
Descrizione - Description - Description.....	3
Etichettatura - Labeling - Etiquetage.....	5
Caratteristiche - Specifications - Spécifications.....	6
Curve di riscaldamento e raffreddamento della cuffia Tube-assembly heating and cooling curves Courbes d'échauffement et refroidissement de la gaine équipée .....	8
Dimensioni - Outline drawing - Schéma dimensionnel.....	9
Fissaggio meccanico / Mechanical mounting / Montage mécanique.....	10
Collegamenti elettrici / Electrical connections / Connexions électriques .....	12
Informazioni sulla sicurezza elettrica - Information about electric safety - Information sur la sécurité électrique.....	15
Informazioni sulla sicurezza termica / Information about thermal safety / Information sur la sécurité thermique .....	16
Modifiche alla filtrazione - Modification of filtration - Modification de la filtration .....	17
Scambiatore di calore - Heat exchanger - Echangeur de chaleur.....	18
Accessori - Accessories - Accessoires .....	23
Simbologia - Symbols - Symboles .....	24

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**Declaration of conformity**

This tube fulfils the essential requirements of the regulation EU 2017/745  
according to standard IEC 60613, IEC 60336, IEC 60522, IEC 60526, IEC 60601-1,  
IEC 60601-1-3, IEC 60601-2-28.

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### Material

Brass body. Filled under vacuum with high dielectric insulating oil, specially processed

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### Compatibility

Suitable for IAE mammography tubes with metal body and grounded cathode

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### Particular configurations

Prearrangement for water circulation with external exchanger.

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### **Intended use**

X-ray units produced by IAE are intended for use in medical diagnostics, in conjunction with equipment having the same purpose.

They are intended to generate X-ray used in radiological equipment with diagnostic functions.

They are not designed to be used in sterile environments, outdoors or in presence of explosive atmospheres.

### **Improper use**

The X-ray tube assembly must not be used if symptoms of mechanical, electrical or radiation-related damage are present.

Manipulation of the security and monitoring devices of the unit is prohibited.

The load limit values listed in this operator manual may not be exceeded.

### **Authorized operators**

In normal operation the user has never access to the tube assembly, because the operation is always controlled by the equipment control board.

For this reason, the tube assembly can only be operated by persons authorized to the use of the radiological equipment where the tube assembly is installed, in accordance with the country specific regulations.

The X-ray tube assembly may be replaced only by personnel trained for service on the specific radiological equipment.

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The windows is equipped with four threaded holes M6 to install the X-ray tube assembly on the equipment and the collimator on the X-ray tube assembly

Minimum number: 4 pcs.  
Minimum strength class: 8.8  
Min engagement length: 14 mm  
Min engagement turns: 14  
Max engagement length: 15 mm  
Max engagement turns: 15  
Max tightening torque: 8.3 Nm

The screws must be secured against possible disengaging by adequate means, according to the expected vibrations during transportation and operation.

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Hazard of loose X-ray tube assembly caused by incorrect tightened screw.  
If you do not obey the above instruction, there is a risk of serious injury.

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Beam limiting devices must interface mechanically with the tube unit so that the continuity of leakage radiation shielding is ensured, and must have a leakage radiation low enough to comply IEC 60601-1-3 or other local regulations, when added to the tube unit leakage radiation.

If a structural element is interposed between the housing and the collimator, the shielding must be integrated with suitable parts to ensure continuity of shielding.

Also, longer screws must be used for mounting, in order to obtain the specified screws engaged length. Compliance with requirements of paragraph Mechanical mounting must be ensured.

A maximum mass of 10 kg is admitted for the collimator.

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Prior than inserting the high voltage cables into the cable sockets, both these and the cable terminals must be thoroughly cleaned. Carbon traces from previous arcings, if any, must be completely removed. Cables must be replaced if they are darkened or swollen near the cable terminals, or if these show any crack.

Apply insulating grease or oil, and additional insulating elements if any, according to the instructions of the cables manufacturer.


It is essential that at the end of the cables installation process, no air or contaminating substance remains in the interface between the mated connectors.

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The standard version of the housing is not provided with watertight protection.

If the housing utilization or its position with respect to the patient or the operator makes it possible spilling of water, disinfectant, organic liquids or other conductive liquids over the housing, the same must be protected with a suitable cover to avoid electrical conduction from the patient or operator to the stator power supply and the thermal safeties.

When this situation occurs, it's absolutely necessary to interrupt the housing use and its power supply, until the liquids have been removed.

The stator power supply must be built so that it is isolated from the mains.

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Two devices are present for thermal protection:

- one bimetallic thermal switch, fitted internally (switches off at 67 °C);
- one bimetallic thermal switch, internally assembled, series connected with stator common cable (switches off at 80°C).

The first thermal switch is accessible to the installer. This one must be connected to the system safety circuits. This device is emergency device. This one must not be intended as a normal operation limiting device. The second thermal switch must be considered as the last safety level; it operates in case of unsuccessful functioning of the other device.

**INSTALLER MUST ALWAYS CONNECT THE FIRST DEVICE EVEN IF THE LAST ONE IS PRESENT !**

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
### Warning

The tube assembly is supplied with an additional 0.03 mm Mo filter available to the installer.

It is the responsibility of the installer to take care that the radiation source assembly complies with the IEC 601.1.3 (CEI EN 60601.1.3).

If needed the supplied additional filter may be used.

*Note: to install the additional filter, please follow the attached instruction.*

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

## Collegamenti - Connections - Connexions

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### Assembling and connection

Fasten the heat exchanger on the equipment by the two M4 tapped holes placed in the box bottom face. For an adequate air flow a 3 cm clearance is required in front of both grids, and the equipment must have air intake and outlet openings, at least 170 square cm wide. Connect the heat exchanger to a 24 VDC  $\pm$  5% , 2.6A power supply by the attached connector, according to the schematics here below.

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
Lay the hoses into the equipment, taking care to avoid any squeezing in all the possible positions of the moving parts.

When assembling is complete, insert the hoses ends on the barbed couplings.

(If the recommended hoses type, Parker Push-lok 801-4, is used, it will only be possible to disconnect the hoses by cutting the ends)

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
## Filling

Prepare a 10% ethylene glycol water solution (quantity depends on hoses diameter and length )

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1. Remove the lateral cover and unscrew the venting plug located in the reservoir base.
2. Unlock and remove the filling plug by pressing the coupling collar downwards. Insert the funnel.
3. Fill the reservoir up to the FILL mark, remove the funnel and reinsert the filling plug.
4. Energize the heat exchanger and monitor the fluid level as it lowers in the reservoir.
5. De-energize the heat exchanger when the reservoir is nearly empty.
6. Remove the filling plug and fill the reservoir again.
7. Repeat the previous steps until the fluid level stabilizes between the MIN and MAX marks.
8. Screw in the venting plug and replace the lateral cover.

*Note: the heat exchanger is operated in a sealed condition (venting plug closed) in order to avoid fluid evaporation, so that a frequent monitoring of the fluid*

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*level is unnecessary. In this configuration, the fluid thermal expansion is compensated by compression of the air volume remaining in the reservoir above the fluid level. This is why the reservoir must not be filled above half of its volume.*

Any excess fluid remaining in the reservoir can be removed by unscrewing the lateral plug in the jumper block and operating the pump shortly.

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Periodically clean the finned array by a vacuum cleaner and a gentle compressed air flow. (Cleaning frequency depends upon the ambient dust concentration). At the same time check the fluid level.

If the flow switch signal is absent, or the tube unit thermal safety devices are activated, all the cooling system elements (pump, fan, finned array, fluid level an hoses) should be checked.

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