Questa documentazione deve essere fornita all'utilizzatore del complesso tubo-guaina  
The contents of this documentation must be transmitted to the user of the tube-assembly  
Le contenu de cette documentation doit être transmis à l'utilisateur de la gaine équipée

Documentazione N° Documentation N° N° de Documentation	Revisione Edition Version	Data di edizione Date of release Date de l'édition	Testo originale Original text Texte original
INFO_TX	D	20.12.12	italiano / italian / italien

## 0. Sommario - Table of contents - Table des matières

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X-ray tubes manufactured by IAE, after a suitable assembling procedure into a tube housing, are intended for use in medical diagnostic, in association with the equipment designed for the same purpose.

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The X-ray tube is a vacuum tube operating from an electrical viewpoint like a diode in conditions of filament saturation.

When high voltage is applied between anode and cathode, the X-ray tube emits X radiation in all the half space above the tangent plane to the anode surface.

For this reason the tube must be used within a shielded envelope, the housing, which limits all the undesired part of X radiation, allowing exit to the useful X-ray beam only, through the housing output port.

The useful beam must be further limited by a fixed or adjustable collimator, which is not a part of the housing and must ensure, together with the housing, the continuity of leakage radiation shielding.

The housing also ensures protection from the high voltage which is applied to anode and cathode by the high voltage cables coming from an elevator transformer with rectifying bridge, usually called generator.

X-ray production requires the use of high energy, most of which is converted to heat dissipated into the tube elements.

The energy sent to the tube initially concentrates on the focal spot (milliseconds), then it distributes over the focal spot track (fractions of a second), then accumulates in the anode mass (seconds to minutes), later is transmitted by heat radiation to the insulating oil which is contained in the tube housing and to the housing metal walls, and from these to the ambient air by natural convection or forced ventilation (fractions of hours).

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These heat exchanges bring the tube elements to extreme temperatures, which must be kept within given limits to avoid destruction of the X-ray tube.

Focal spot and focal spot track temperatures are controlled by power (kV x mA) and exposure time (sec), which are supplied by the generator under the limitation of maximum load factors, as calculated from the rating charts.

Anode temperature is governed by a correct exposure to cooling delays ratio, as a function of the exposure energy.

The control of cooling delays must be ensured by the equipment resident software, by accumulated heat simulation circuitry or, in absence of these, by a careful and conservative programming of working schedules made by the user, on the basis of anode heating and cooling curves.

Housing temperature is also controlled by work to cooling delays ratio, by the same means described for the anode. In this case the duty cycles are calculated over longer intervals, of the order of half a working day, on the basis of housing heating and cooling curves. Energy computation must include the energy which is sent to the stator for rotor acceleration, which can be of the order of 10 kJ per cycle in high speed operation.

Up to three safety devices are present in the tube unit for the housing temperature only: an external thermal switch, an internal one not accessible to the user and a microswitch on oil expansion bellows.

These devices are not to be considered as a substitute to automatic or manual programming of cooling delays, but as last safety levels never to be reached in normal operation.

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**possible safety device exist in the tube unit for the other limiting temperatures, concerning focal spot, focal spot track and anode. For these, patient safety and tube integrity must be ensured by a correct control of energy cycles from the generator.**

Anode current as a function of filament current is given by the emission curves.

For construction reasons the tolerance on the filament current value which gives a determined value of anode current is comparatively wide (approx. 3% of filament current). For this reason the emission curves can be used as an initial operation set up, but cannot replace an individual tube calibration at each installation.

Pre-programmed systems containing information on filament emission must also allow for the same range of filament calibration.

The useful beam is filtered by the tube unit structural materials (glass, oil, plastic materials). However this inherent filtration is not sufficient to limit the low energy fraction of X radiation, the so called soft X-rays, which are harmful to the patient and do not contribute to image creation. This happens even more with Beryllium window tubes, typically mammography tubes, which emit a high quantity of soft X-rays.

For this reason an additional filtration must be placed on the useful beam, ensuring the compliance to the pertinent regulations while keeping account of all the interposed materials (collimator, compressor etc.)

A part of added filtration is normally present on the tube unit, but the obtaining of a sufficient total filtration is a task of the installer.

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tient. Since the diagnostic benefit to harm ratio can only be evaluated by a physician, X-ray equipment are to be operated exclusively by competent personnel under the guide of a radiologist.

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The X-ray tube operates in conjunction with the following components of the radiological equipment, whose compatibility must be checked:

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As a prerequisite to use, the housing must be safe, according to the pertinent regulations, for what concerns leakage radiation protection, tube unit explosion and mechanical safety. When selecting an X-ray tube to replace an original tube in an existing housing, the following elements must be taken into consideration:

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The distance of the glass bulb central part from housing metal walls should be preferably between 12 and 15 mm, and in any case not lower than the original tube distance.

The distance between the glass bulb conical surface behind the target and the stator insulating sleeve must be not less than 2 mm, to allow for thermal expansion between tube and housing without danger of glass stress.

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#### 4.1.2 Supporto meccanico - Mechanical support - Support mécanique

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The tube anode support must give a reliable and rigid grip on tube anode end.

If the housing internal structures are suitable, the tube should be held by the anode end only, to avoid dangerous stresses on the glass bulb because of misalignments and thermal

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expansions.

Information on specific housing types, as well as adaptation pieces, can be obtained from IAE Spa.

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Standard rotor tubes are suitable for operation with stators having a stator core inner diameter not larger than 62 mm.

In case of wider stators, tubes with larger rotor diameters must be selected.

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The oil gap between the window port and the glass bulb should be between 3 and 5 mm. A lower gap gives an insufficient insulation, a higher one increases the tube assembly inherent filtration.

In the latter case, a higher oil gap may be accepted, provided that the increased oil thickness is compensated by a suitable variation of additional filtration.

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If the tube focal spots position with respect to the housing is not the same as the original tube, the effects of this difference on film coverage and collimator field limitation must be taken into consideration.

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The insulating oil dielectric strength must comply with the pertinent regulations.

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Both high frequency and mains frequency generators can be used. The X - ray tube is a passive component not being able to limit the energy sent by the generator.

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For this reason, the generator - control board system must comply to the pertinent regulations, requiring that in any failure condition an uncontrolled energy supply to the tube is prevented as this eventuality could cause the tube unit explosion with severe danger to the patient and the user.

When designing a new equipment or replacing an existing tube with a new tube of different make, it must be checked that the filament heating circuit can supply the filament voltage and current as required by filament characteristics and emission curves, including the calibration range required by filament current tolerances indicated on the same curves.

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### 5. Trasporto - Transportation - Transport

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X ray tubes, both new and returned, must always be shipped in the original packing, after checking that the packing is in good conditions.

While the housing structure is made by sturdy metal parts, the fragile part is the glass insert inside.

It is thus important that during all the transportation phases the tube be not subject to shocks and be kept in vertical position, the situation in which the glass part is less sensitive to shocks.

Should the packing be open for inspection and then undergo further transportation, care must be taken that all the packing parts are restored in their original position.

Upon receipt of the tube, the shock and tilt indicators should be checked for activation, and the subsequent actions taken according to the indicators instructions.

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Tubes are to be stored in a non-oxidizing, non-corrosive, dry environment, within the temperature limits indica-

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ted in the data sheets.  
Should the tubes be kept in stock for more than a year, the factory or the local representative should be contacted in order to repeat the tube test.

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Reloading a tube assembly with a new tube must be done by qualified personnel, according to specific instruction which can be obtained from the manufacturer.

This operation requires handling of the glass bulb, which, as a vacuum tube, can implode in case of shocks or damages.

The operator must thus protect himself by a face shield and heavy gloves from tube unpacking to tube insertion into the housing.

As general information, the following issues must be taken into consideration when reloading a tube, in addition to the issues listed in the *Interfacing and compatibility* section.

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All the parts contributing to leakage radiation shielding must be checked for presence and compatibility. Damaged parts must be replaced, or the housing must be discarded.

In case of doubt concerning shielding efficiency, a leakage radiation survey must be performed before the installation.

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All the insulating parts (high voltage cable sockets, stator sleeve, anode insulator) must be checked for absence of surface arcing or high voltage puncturing.

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The correct values of resistance and insulation to ground must be present at the stator terminations.

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The housing internal surfaces must be clean and free from carbon residuals or old oil, and from sharp corners which could start high voltage arcs.

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---

**7.5 Supporti meccanici - Mechanical supports - Supports mécaniques**

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The mechanical supports must be compatible with the tube anode termination and give an efficient tightening.

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The insulation distances between the parts connected to high voltages and the housing grounded parts must not be smaller than the original setting. The distance between the glass bulb central part and the housing metal walls should be preferably between 12 and 15 mm.

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Check that the gap between the glass bulb, in the area behind the target, and the stator sleeve is enough to allow tube and housing thermal expansions without stressing the glass.

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In case that because of the housing internal structure the tube needs to be supported also by the cathode end, check for misalignments or forces that may apply dangerous stresses to the glass.

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For a proper insulation the minimum oil gap between glass and window is 3 mm. Significant variation of oil gap referred to the original assembly must be compensated by a variation of added

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In order to centre correctly the focal spots on centre of window port, the following alignments must be obtained on the tube:

- cathode centred with respect to the windows axis parallel to the housing main axis.
- external edge of the target conical surface 2 mm off centre towards housing anode end.

A collimator must be used for a correct alignment.

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Care must be taken in avoiding stresses on the filaments feed throughs, which could cause glass cracks and vacuum leaks.

Check that the filaments are connected correctly to ensure a correct selection of focal spots in operation.

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The anode and filaments connections must follow the original paths in terms of distances from grounded parts, bearing in mind that these connections will assume respectively the positive and negative high voltages.

Avoid longer than needed connections, which could approach grounded parts by electrostatic attraction when applying high voltage to the tube unit.

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Oil filling must be performed under vacuum in order to eliminate air residuals, and preceded by a drying process to remove moisture from the housing internal surfaces.

Oil dielectric strength must comply to standards for processed insulating oil,

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and be checked on an oil sample taken from the housing after filling, to take account of possible oil contamination by impurities or moisture residuals from the housing surfaces.

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The tube assembly must undergo a high voltage test before installation, to check the tube stability and the housing insulation efficiency.

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Check by the measurement of the filament transformer primary current and coil ratio, that the filament preheating current is between 2 A and 3 A. Too high a preheating current makes the filament evaporate, metallizes the glass and shortens the tube life.

Check that in radiography preparation the large and small filaments are selected correctly, to prevent focal track overloading in case of wrong focal spot selection.

**WARNING** to perform this check the high voltage to the tube must be inhibited, to avoid the risk of being irradiated from a short distance.

Make a shortened rotation boost to check the correct rotation direction, i.e. anode rotating CCW when looking at the anode from the cathode side.

Check that the rotor full speed is attained by the normal boost cycle. In RTM tubes the rotor speed can be measured by the frequency of the rotor feedback voltage in the stator windings.

To perform the measurement, the starter has to be disconnected from the stator, allowing the rotor to rotate by inertia. In these conditions a voltage of approx. 400 mV, at the rotation frequency, is present between the stator common wire and each of the other two wires. ( $RPM = 60 \times f$  (Hz)). In the other tubes the speed can be determined by a vibrometer, a stro-

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boscope or a reed frequency meter.

A non quantitative indication can be obtained by the rotation noise. If the boost cycle does not attain full rotor speed, the noise frequency will increase during rotation stand-by.

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If the last tube test was earlier than two months, perform the seasoning procedure on the installation according to the supplied instructions.

The same procedure must be performed if the tube has a tendency to make arcs.

In any case perform the warm up section of this procedure.

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Perform the calibration according to the equipment instructions.

To prevent tube unit overheating and to avoid and excessive lengthening of the calibration procedure, it may be useful to switch the stator power supply to the stand-by voltage.

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In normal operation the user has never access to the tube, because the tube operation is always controlled by the equipment control board.

For this reason, specific instructions for the tube cannot be given. Reference should be made to the equipment instructions.

Below some general purpose indications are given, being relevant or safety of use or tube life.

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Applying high loads to the tube with cold anode can prematurely wear the focal track, because the anode material has a fragile behaviour at room temperature, and thus its surface is more easily cracked under the thermal cycles induced by the load application.

While feasibility or emergency reasons are often predominant in diagnostic use, it is however important for tube life to perform a warm-up, as it is usual in computed tomography installations, if the tube has been left inoperative for more than 3 hours and the tube is operated at the highest powers. If the warm-up procedure is not included in the equipment software, the necessary data can be obtained by the first part of tube seasoning procedure here included.

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While the control board normally limits the peak radiographic power (maximum loads) this not always happens for the anode and housing thermal limits (accumulated energy and continuous power dissipation). However these are extremely important for a safe use of the equipment.

In absence of an equipment resident software having this purpose, the daily working schedule must be determined, with the help of the installation engineer if needed, so that:

- the energy of one examination is lower than the maximum anode heat content
- the average power of consecutive exams (work + delays) is lower than the maximum anode dissipation
- the average power of half a working day is lower than the maximum housing dissipation

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If the equipment has remained inoperative:

- for more than a month, perform the warm-up cycle
- for more than 6 months, perform the seasoning procedure.

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Normally tube failures (high voltage arcs, filament cut off, stator failures) are detected by the equipment safety devices.

However it is important to act at a prevention level every time any abnormal noise coming from the tube unit or other parts of the high voltage circuit (cables, generator) may give suspect of a state of failure.

Under these circumstances **suspend exposure immediately, cut off power to the equipment and remove the patient**, then contact the service engineer and follow his instruction.

Also perform an immediate action, either directly or through the service engineer, if any suspended mass (housing or collimator) shows an unstable fastening.

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In the occasion of ordinary equipment maintenance, the following checks are required by the tube unit:

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By means of an oscilloscope and a high voltage probe, either invasive or not, that no part of the high voltage circuit, including the tube, show a tendency to arcs.

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The tube contains the following materials, whose disposal must be done in respect of the local regulations:

- Tungsten, molybdenum, rhenium, nickel, copper, silver;
  - Glass;
  - Iron, aluminium and copper alloys;
- Spent tubes can be sent to the factory for disposal

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In addition to the above, mammography tubes contain an output window made by beryllium, which is highly toxic if ingested.

This part must be in no way machined or worked in a manner suitable of producing powders or fumes, sand blasted or acid etched.

We recommend to send the spent mammography tubes to the factory, leaving them untouched.

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IAE performs all the possible actions, at design, manufacturing and test level, to make safe the use of the tubes of its production.

However there exist residual risks connected to the physical principles of operation, which are characterized by high voltages, temperatures and energies being present in the x-ray tube, and by a limited and random life of the components.

These risks are to be evaluated vs. the diagnostic benefits and be well known to those who must take precautions on the subject.

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---

### 12.1 Radiazioni - Radiations - Radiations

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X-rays emitted by the tube vary in intensity and energy according to the equipment setting. While being useful from a diagnostic

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viewpoint, X-rays are anyway harmful to the patient.

*The equipment use must thus be strictly limited to qualified personnel under the guide of a radiologist, who can determine the correct exposure parameters and establish a risk to benefit ratio.*

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The glass vacuum tube can be broken because of shocks or manufacturing defects, causing the bulb implosion and glass fragment projection against the reloading operator or any person handling the glass tube.

*Anyone handling the glass tube must use protective face shield and gloves.*

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In case of alteration of the focal spot size compared with those required by the pertinent regulations, the image quality can be impaired, causing possible diagnostic errors.

*The user must perform periodical image quality tests, as required by the pertinent regulations.*

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All the seasoning procedure should be performed by the large focal spot. Stator boost cycles must be inhibited during serial exposures in order to prevent tube unit overheating and stator damage.

Perform the following steps by selecting the technique as shown in the tables.

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Perform the required number of exposures with kV, mA, exposure time and cooling delays as shown in the table.

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Starting from the listed initial kV values, make a set of exposures with the required kV, mA, exposure time and cooling delays. At the end of the series increase kV by the required step and repeat the series. When the final kV value is attained, make the final number of exposures as shown in the table.

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Starting from the initial kV value, perform fluoroscopy at the listed kV, mA, fluoroscopy time and cooling delay. Step up kv by the listed value and repeat fluoroscopy. When the final kV value is attained, make the final number of fluoroscopy cycles.

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Starting from the initial kV value, make radiographic exposures at the listed kV, mA, exposure time and cooling delay. At the end of the series increase kV by the required step and repeat the series. When the final kV value is attained, make the final number of exposures as shown in the table.

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If normal operation does not attain the maximum kV values, the seasoning cycles can be limited to the kV values of actual use.  
The warm-up phase is also intended as a daily warm-up.

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